

SurfaSoft[®]

Monofilament, non-adherent, micro-perforated, protective contact layer

- Fixation dressing for skin grafts / dermal matrices (i.e. MatriDerm[®])
- Fixation dressing for RECELL[®] Spray-On Skin
- Interface dressing for NPWT
- Wound contact layer



SurfaSoft®

Monofilament, non-adherent, micro-perforated, protective contact layer

Product Description

SurfaSoft® is a unique transparent dressing made of monofilament polyamide. The rolled threads create a mesh with a pore size of 90 microns, producing a strong, smooth and transparent material with an open draining surface of 41%.

SurfaSoft® can be used as fixation material for skin grafts and/or dermal matrices, as an interface for negative pressure wound therapy (NPWT) and as a wound contact layer.

This unique porous pattern allows exudate to evaporate and drain without drying the wound surface, while topical preparations used to treat burns (post-surgery) can pass through to the skin graft.

'State-of-the-art' transparent fixation dressing for skin grafts (sheet, mesh or meek) and dermal matrices i.e. MatriDerm®

SurfaSoft® is a strong, flexible, inert, non-adherent dressing that reduces the friction and shear forces when used for fixation of skin grafts.

SurfaSoft® ensures a tight skin graft and wound fit for optimal graft ingrowth. The SurfaSoft® dressing fixes skin grafts and promotes wound drainage, while minimising bacterial adhesion. When used with SurfaSoft®, absorbent wound dressings can be removed painlessly and as often as needed. SurfaSoft® allows for easy wound inspection, can remain in situ for several days (or until wound healing) and can be removed without damaging the underlying tissue. With SurfaSoft® dressings, skin grafts in difficult areas, e.g. the armpit, neck, head or back, have a take rate as high as 100%.



SurfaSoft® application and fixation for the protection of autografts.



SurfaSoft® day 2 post grafting and remains in place for 5 to 6 days



Removal of SurfaSoft® on day 6 post autografting, with 100% graft take observed

SurfaSoft®

Monofilament, non-adherent, micro-perforated, protective contact layer

The ideal interface when applying negative pressure wound therapy

SurfaSoft® can be used as a protective contact layer for large soft tissue defects, large abdominal defects, protecting fragile and vital structures such as the gastrointestinal tract, heart, exposed tendons, ligaments, prosthetics and more.

SurfaSoft® also limits granulation tissue ingrowth in polyurethane foam and makes dressing changes easier and less painful for the patient.

The perfect wound contact layer

As a non-adherent interface SurfaSoft® can be used to protect vulnerable wound surfaces and promotes undisturbed wound healing by preventing adhesion to the wound bed (and granulation tissue damage) or to the newly formed epithelium. SurfaSoft® maintains an optimal moist environment for wound healing. SurfaSoft® can be used safely with most absorbent dressings thanks to the inert contact material. Dressing changes are also much easier to perform and are less painful for the patient.

Product features and benefits

Easy use

- Easy to cut to size
- Does not stick to gloves
- Very easy to apply and to model in cavities and underminings
- Easy to remove without damaging the newly formed epithelium or granulation tissue

Tissue-friendly

- Non-adherent
- Non-irritating
- Hypoallergenic
- Latex free



Product applications

- Autografts and allografts
- Donor sites
- NPWT indications
- Dermal Matrices
i.e. MatriDerm® Dermal Matrix
- RECELL® Spray-On-Skin
- Partial & full-thickness wounds
- Traumatic and chronic wounds
- Clean closed surgical incisions
- Dermatologic lesions

SurfaSoft®

Monofilament, non-adherent, micro-perforated, protective contact layer



DESCRIPTION

SurfaSoft® 10 cm x 15 cm

SurfaSoft® 20 cm x 33 cm

ORDERING CODE

1015S

2033S

PACKAGING

Box, 10 pcs

Box, 25 pcs



HAROMED B.V.

Beukenlaan 21 | 9051 Gent | Belgium | Tel +32 9 326 05 10
info@haromed.com | www.haromed.com



Distributed by **Revolution Surgical Pty Ltd**
4/22 Rowood Road, Prospect NSW 2148
Email: info@revolution-surgical.com
Website: revolution-surgical.com